

Using 7T Multimodal Data

Junqian (Gordon) Xu, An (Joseph) Vu

HCP OT 1-3: UMN (CMRR), Oxford, and WU

<https://dl.dropboxusercontent.com/u/23744124/HCP7T.pdf>

Learning objectives

- Advantages and challenges of 7T neuroimaging
- MR techniques developed at CMRR, UMN to meet the challenges
- BOLD contrast mechanism and its field dependence
- HCP 7T protocol
- HCP 7T data preliminary analysis results

<https://dl.dropboxusercontent.com/u/23744124/HCP7T.pdf>



Siemens 7T AS 90/60 cm bore
SC72CD gradients

$G_{\max} = 70$ mT/m (for diffusion)

$G_{\max} = 42$ mT/m (for imaging)

Slew rate = 200 T/m/s (for imaging)

32ch Nova head coil

- 3T Phase 2: Aug 2012 - 2015

1,200 participants, 300 twin pairs

- 7T Pilot: Aug 2012 - late 2013
- 7T Phase 2: late 2013 – 2015

a subset of 200 HCP subjects
scanned at CMRR, UMN

Advantages of 7T MRI

- **Increased SNR traded for spatial resolution**
- **Increased intrinsic BOLD sensitivity**
- **Increased BOLD specificity (in certain regions)**

MR sensitivity ~ field strength

$$S/N \sim B_0^{3/2}$$

system noise
dominant

SNR gain $_{7T/3T}$: 3.75!

7T (297 MHz)

3T (123 MHz)

$$S/N \sim B_0$$


sample noise
dominant

SNR gain $_{7T/3T}$: 2.4

Challenges of 7T MRI



- **Relaxation: longer T1 and shorter T2/T2***
- **Transmit: inhomogeneous B_1^+ field**
peak power limit
specific absorption rate (SAR)
- **Main field: inhomogeneous B_0**
- **Motion: increased sensitivity (bulk or physiological)**

Importance of Improving B_1^+ Homogeneity at 7T



- Over flipping in the center of the brain
- Low flip angle at lower part of the brain (temporal lobe, cerebellum, brain stem)

Dielectric material to improve B_1^+





Dielectric pads
(CaTiO₃)

Over flipping in the center of the brain remains


Time and Phase shifted MB pulses to reduce peak power and specific absorption rate (SAR)

Peak Power and BW with 25% temporal shift between bands




Motion during the reference/calibration scans (higher resolutions + MB+IPAT)

standard EPI Reference, 5 different scans



GRE Reference acquired once at the beginning

GRE reference





Advances in HCP 7T MRI acquisition


128 slice limit eliminated

Dielectric pads (CaTiO₃)

Motion robust GRE reference scan



Low SAR fat suppression




Learning objectives

- Advantages and challenges of 7T neuroimaging
- MR techniques developed at CMRR, UMN to meet the challenges
- BOLD contrast mechanism and its field dependence
- HCP 7T protocol
- HCP 7T data preliminary analysis results

Questions?

BOLD contrast mechanism

- deoxy-Hb (paramagnetic)
- oxygen metabolism
- cerebral blood flow (CBF)
- cerebral blood volume (CBV)




Barth and Poser, 2011

BOLD signal contribution

- Intravascular (IV)
- Extravascular (EV)

Arrangement of pial veins





Auer and Loew, 1983






capillary ~ 5 μm
venule \updownarrow
vein ~ 500 μm - mm

Extravascular (EV) contribution





Superlinear field dependence of BOLD sensitivity



- Extravascular (EV) contributions all increase as field increases
- % intravascular (IV) contribution (esp. large vein) decrease as field increases
- % contribution from microvasculature increases as field increases







Reduced partial volume increases BOLD sensitivity



Resting State Networks at 7T vs 3T Connectome

7 Tesla, 1.25 mm Isotropic Res.
MB=3, IPAT=3, TR=2 s

3 Tesla, 2 mm Isotropic Res.
MB=8, IPAT=0, TR=0.7 s



Learning objectives

- Advantages and challenges of 7T neuroimaging
- MR techniques developed at CMRR, UMN to meet the challenges
- BOLD contrast mechanism and its field dependence
- HCP 7T protocol
- HCP 7T data preliminary analysis results

Questions?

Criteria for HCP 7T Acquisition

- whole brain coverage
- acquisition & image reconstruction stability
- relative scalability (200 subjects)
- spatial resolution
- temporal resolution (fMRI)

HCP 7T session structure

- **Session 1 (15 min rfMRI, 30 min movie)**
- **Session 2 (15 min rfMRI, 40 min dMRI)**
- **Session 3 (15 min rfMRI, 30 min retinotopy)**
- **Session 4 (15 min rfMRI, 30 min movie)**

fMRI (total): 1hr resting, 1hr movie, 0.5hr retinotopy

No anatomical MRI at 7T

HCP fMRI Protocols

HCP 3T

2 mm iso

72 slices / MB8

iPAT = 0

PE = LR/RL

TE/TR = 33/720 ms

HCP 7T

1.6 mm iso


85 slices / MB5

iPAT = 2, 7/8 partial Fourier



PE = AP/PA

TE/TR=22/1000 ms

Number of ICA components *reduced* at higher spatial resolutions



Cerebral Cortical Thickness vs Imaging Resolution





- 3T: 2.0mm resolution, 1 frame / 0.72s
- 7T: 1.6mm resolution, 1 frame / 1.0s
 - 7T fMRI goal was to somewhat increase spatial resolution while gaining a substantial amount of CNR

Resting State Denoising


- Very gentle highpass filter removes largely linear trends
 - sigma of 1000s, runs are 900s in length
- ICA with up to 250 components
 - automatically determined (in general higher SNR or more structured signals will lead to higher dimensionality)
- FIX classifies ICA components into signal and noise categories
 - In this case, because FIX was being retrained for 7T data, all components were hand classified
- 24 motion parameters and noise component timeseries are regressed out of the data

3T vs 7T: rfMRI CNR







Number of Total ICA Components




Number of Noise ICA Components




Number of BOLD Signal ICA Components



3T vs 7T: Dense Connectivity



3T vs 7T: Dense Connectivity



Learning objectives

- Advantages and challenges of 7T neuroimaging
- MR techniques developed at CMRR, UMN to meet the challenges
- BOLD contrast mechanism and its field dependence
- HCP 7T protocol
- HCP 7T data preliminary analysis results

Questions?

HCP dMRI Protocols

HCP 3T

1.25 mm iso

111 slices / **MB3**

iPAT = 0, 6/8 partial
Fourier

PE = LR/RL

TE/TR=89/5500 ms

b = 1, 2, 3k HARDI

HCP 7T

1.05 mm iso

132 slices / **MB2**

iPAT = **3**, 6/8 partial
Fourier

PE = AP/PA

TE/TR=71/7000 ms

b = 1 & 2k HARDI

7T improved spatial resolution with less blurring


7T
iPAT3, MB 2, 70 mT/m

3T HCP Skyra
MB 3, 100 mT/m

b=0 s/mm²

b=1000 s/mm²



b=2000 s/mm²




DTI Maps at 7T vs 3T Connectome

7T, 1.05 mm iso, MB2, IPAT3, 143 dirs
TR = 7000ms, PF 6/8, TA: 38mins
b = 1 & 2k s/mm², G = 70 mT/m
time shift = 1920 μ s, dielectric pads




3T, 1.25mm iso, MB3, IPAT1, 288 dirs
TR = 5500ms, PF 6/8, TA: 60mins
b = 1k, 2k, & 3k s/mm², G = 100 mT/m






7T HCP Fractional Anisotropy Map



Superficial dense white matter fibers (Postmortem human brain)



3T vs 7T HCP dMRI Bedpost







Can we analyze the datasets jointly & retain features from both?

(Sotiropoulos et al ISMRM # 562)

HCP 7T and 3T Diffusion

(Sotiropoulos, Vu et al)




Learning objectives

- Advantages and challenges of 7T neuroimaging
- MR techniques developed at CMRR, UMN to meet the challenges
- BOLD contrast mechanism and its field dependence
- HCP 7T protocol
- HCP 7T data preliminary analysis results

Questions?

fMRI tasks 3T


- Working memory
- Gambling
- Motor
- Language
- Social Cognition
- Relational Processing
- Emotion Processing




Barch et al. *Neuroimage* 2013

fMRI tasks 7T

- Retinotopy



- Natural movies w/ sound



Movie Clips 4 x 15 min runs

Hollywood (Cutting 2012)

- Inception Run 1
- The Social Network
- Ocean's Eleven



- Home Alone Run 3
- Erin Brockovich
- Star Wars Episode V

Creative Commons (vimeo.com)



- Two Men Run 2
- Welcome to Bridgeville
- Pockets – Short Film
- Inside the Human Body

- Off the Shelf Run 4
- 1212
- Mrs. Meyer's Clean Day
- Northwest Passage (montage)

Movies elicit broad and robust brain activations



4 X 15 mins
2 mm iso
iPAT2 MB5




Decoding the neural coding

Movie stimuli:

- 14 movie clips
- 55 min total

First 3 Semantic
Principal Components in

- **purple** (buildings and vehicles),
- **green** (faces),
- **yellow** (animals).



Learning objectives

- Challenges and advantages of 7T Imaging
- BOLD contrast mechanism and its field dependence
- MR techniques developed during HCP 7T at CMRR
- HCP 7T protocol
- HCP 7T preliminary analysis results


Questions?

Retinotopy Introduction

- Visual cortex can be divided into distinct areas based on mapping out the representation of the visual field
- Use fMRI to determine retinotopy; use 7T for higher SNR

- What is the goal?

- Find the visual areas (so you can use them as ROIs?)
- Discover new ones?
- Use the large subject pool ($n = 200$) to get population statistics?




Larsson 2006

- A little history / background:

- First set of pilots played with experimental design
- Second set of pilots settled on a particular design
- Experiment is run using MATLAB + Psychtoolbox on Mac OS X


Retinotopy Experiment

- Stimulus type
 - “mashfast” (multiscale objects on pink-noise background, 15-Hz)
 - Evokes larger responses, improves reliability of phase estimates
- Aperture and run design (each run is 300 s)
 1. CCW. 22 s; 8 cycles * 32 s; 22 s
 2. CW. "
 3. EXPAND. 22 s; 8 cycles * (28 s + 4 s); 22 s
 4. CONTRACT. "
 5. BARS. 16 s;
 - R, U, L, D (each 28 s + 4 s); 12 s;
 - UR, UL, LL, LR (each 28 s + 4 s); 16 s
 6. BARS. "
- Task: detect change in color of fixation dot (every 1–5 s)
- Eye-tracking
- Field-of-view: $\sim 15^\circ$ stimulus diameter
- fMRI: 1.6 mm, 1 s




Angle

Wedges & Rings – attend to stimuli




Angle

Right & Up Bars – attend to stimuli




Eccentricity

Wedges & Rings – attend to stimuli



Eccentricity

Bars – attend to stimuli



Learning objectives

- Advantages and challenges of 7T neuroimaging
- MR techniques developed at CMRR, UMN to meet the challenges
- BOLD contrast mechanism and its field dependence
- HCP 7T protocol
- HCP 7T data preliminary analysis results

Questions?

HCP Minimal Preprocessing Pipeline for 7T data


- Pipeline code: upcoming 3.6.0 release
- Requires 3T structural (T1w/T2w)
- Siemens SC72CD unwarp file (gradient non-linearity correction)
- ICA-FIX retrained for 7T (fMRI)
- Latest FSL/eddy version (dMRI)

**More details in the upcoming HCP 7T data
release documents**



Acknowledgments

An (Joseph) Vu



Essa Yacoub



Kendrick Kay



Mathew Glasser



Tim Behrens



Stam Sotiropoulos



Saad Jbabdi



Jesper Andersson

Center for Magnetic Resonance Research UMi



Kamil Uğurbil



Steve Smith



Karla Miller



Eddie Auerbach



Steen Moeller



Christophe Lenglet